



A VERSATILE, BROAD SPECTRUM DISINFECTANT FOR USE IN CRITICAL CLEANROOM AREAS

Klercide™ 70|30 IPA

70% v/v Isopropyl alcohol blended with Water for Injection

Klercide 70|30 IPA blended with Water for Injection is available in a variety of sizes and formats for optimum performance depending on application method required.

Klercide 70|30 IPA is:

- ▲ 70% v/v Isopropyl Alcohol blended with Water for Injection
- ▲ Available in sterile and non-sterile, filtered options
- ▲ Available double or triple bagged
- ▲ Supplied in our patented and validated SteriShield Delivery System

The SteriShield Delivery System (SDS) features:

- ▲ A closed trigger system providing a completely sealed bottle
- ▲ Protection for the integrity of the contents throughout use
- ▲ Validated in-use shelf life of three months
- ▲ A fully adjustable spray pattern to ensure optimum surface coverage

PRODUCT QUALITY ASSURANCE

Manufactured in our state-of-the-art facilities to the principles of GMP, Klercide 70|30 IPA is safe for use in critical aseptic core manufacturing areas.



Klercide™ 70|30 IPA



70% v/v Isopropyl alcohol blended with Water for Injection

Why use Klercide 70|30 IPA?

OPERATOR SAFETY

Klercide 70|30 IPA is available in Ecolab's proprietary SteriShield Delivery System (SDS).

Compared to an aerosol spray, the larger droplet size from the trigger spray reduces the level of atmospheric alcohol created at the time of use, consequently reducing both the risk of excessive vapor inhalation by operators and explosion risk from aerosols which are packaged under pressure.

REDUCED COST OF USE

Klercide 70|30 IPA is available in Ecolab's proprietary SDS, as it provides the same "closed system" benefits of an aerosol product, but may reduce costs by eliminating expensive, pressurized metal containers and their associated waste disposal requirements, as well as waste of unused product that cannot be dispensed.

From purchase to disposal, the SDS can offer a significant advantage by reducing the total cost of use.

DESIGNED FOR CLEANROOMS

Klercide 70|30 IPA is manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klercide 70|30 IPA is manufactured to the highest quality and are compliant for use in critical areas.

Klercide 70|30 IPA is:

- ▲ Made with 70% v/v Isopropyl alcohol blended with Water for Injection
- ▲ Endotoxin tested to an approved pharmacopoeia method
- ▲ Passes AOAC disinfectant standard for bactericidal and fungicidal efficacy
- ▲ Filled and capped in an ISO Class 7 cleanroom
- ▲ Supported with Certificates of Analysis and Sterility with every batch
- ▲ Supported with full documentation

MICROBIOLOGICAL TEST RESULTS

EFFICACY	CONTACT TIME
Disinfection	5 min
Virucidal	5 min
Fungicidal	10 min and 5 min*
Non-Food Contact Sanitizer	3 min

*Please see label for specific claim

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences
www.ecolab.com/pharmaceutical

© 2019 Ecolab USA Inc. All rights reserved. 17JULY19/NA

Safety. Efficiency. Quality.

APPLICATION AREA

For use on floors, walls and hard surfaces.

INSTRUCTIONS FOR USE

Spray: For optimum results hold bottle upright and spray from a distance of 6" to 8". Always close the nozzle after use.

Wetting: For optimum results apply solution to wet all surfaces thoroughly. All surfaces must remain wet for the required contact time indicated in the directions for use on the product label.

INGREDIENTS

70% v/v Isopropyl alcohol blended with Water for Injection.

REGISTRATION DETAILS

EPA Reg. No. 1677-249
Canadian DIN 02447789

PACKAGING INFORMATION

Klercide 70|30 IPA blended with WFI

UNIT SIZE	UNIT OF SALE
1 L spray	6 x 1 L
16 oz spray	12 x 16 oz
32 oz spray	6 x 32 oz
1 Gallon capped	4 x 1 gal.

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

SHELF LIFE

Unopened, this product has a three year shelf life.



Always supplied with the SDS - the patented and validated trigger spray system

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.