



EFFECTIVE CLEANING SOLUTION FORMULATED SPECIFICALLY FOR CLEANROOM USE



Klercide™ Neutral Detergent

A low foaming, non-ionic detergent solution

An ideal solution that can be used on cleanroom surfaces and production areas. For use in general maintenance and clean-up, as well as residue removal.

Klercide Neutral Detergent is available in multiple formats, including a one gallon capped container and ready-to-use 32 oz. spray bottle.

Klercide Neutral Detergent is:

- ▲ Effective against oily or sticky spillages
- ▲ Designed to clean in critical areas
- ▲ A ready-to-use formulation
- ▲ 32 oz. spray bottle supplied in our patented and validated SteriShield Delivery System

The SteriShield Delivery System (SDS) features:

- ▲ A closed trigger system providing a completely sealed system
- ▲ Protection for the integrity of the contents throughout use
- ▲ Validated in-use shelf life of three months for 32 oz. spray bottle
- ▲ A fully adjustable spray pattern to ensure optimum surface coverage

PRODUCT QUALITY ASSURANCE

Manufactured in our state-of-the-art facilities to the principles of GMP, Klercide Neutral Detergent is safe for use in critical aseptic core manufacturing areas.



Klercide™ Neutral Detergent



A low foaming, non-ionic detergent solution

Why use Klercide Neutral Detergent?

UNIQUE FORMULATION

Klercide Neutral Detergent is a low foaming, non-ionic detergent with a neutral pH. It is a mix of fatty alcohol ethoxylate that is produced, 0.2 micron filtered, filled and capped in an ISO Class 7 cleanroom.

Non-corrosive and highly-effective when used directly on most non-porous surfaces leaving minimal residues for easy removal.

SIMPLE

Ready-to-use format reduces handling and need for dilution. No filtration required. This results in cost effective, simple solution for a wide range of cleanroom contamination control challenges.

DESIGNED FOR CLEANROOMS

Klercide Neutral Detergent is manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klercide Neutral Detergent is manufactured to the highest quality and are compliant for use in your critical areas.

Klercide Neutral Detergent is:

- ▲ Double bagged
- ▲ Supported with Certificates of Analysis and Sterility with every batch



Klercide Neutral Detergent
32 oz. SDS spray



Klercide Neutral Detergent
1 gallon capped

Effective. Cleansing. Direct.

APPLICATION AREA

For use on floors, walls and hard surfaces.

INSTRUCTIONS FOR USE

Spray: For optimum results hold bottle upright and spray from a distance of 6 to 8 inches. Always close the nozzle after use. Dispose of packaging only if emptied and closed.

Wetting: For optimum results ensure complete wetting of the surfaces using the spray and wipe technique.

INGREDIENTS

Fatty alcohol ethoxylate in Water for Injection.

PACKAGING INFORMATION

Klercide Neutral Detergent

UNIT SIZE	UNIT OF SALE
32 oz spray	6 x 32 oz
1 Gallon capped	4 x 1 gal.

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

SHELF LIFE

Unopened, this product has a two year shelf life.

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences
www.ecolab.com/pharmaceutical