



**EFFECTIVE CLEANING SOLUTION FORMULATED SPECIFICALLY FOR CLEANROOM USE**



# Klercide™ Neutral Detergent Unit Dose Concentrate

A low foaming, non-ionic detergent solution

**An ideal solution that can be used on cleanroom surfaces and production areas.**

Klercide Neutral Detergent is available in a variety of sizes and formats for optimum performance, including a soluble unit dose concentrate format that makes up to 1 gallon of ready-to-use product. It is ideal for use in large areas such as walls, ceilings and floors.

Klercide Neutral Detergent Unit Dose Concentrate:

- ▲ Easily creates the correct, validated concentration reducing the likelihood of a dilution error
- ▲ Requires minimal handling to dilute
- ▲ Effective against oily or sticky spillages
- ▲ Designed to clean in critical areas

**PRODUCT QUALITY ASSURANCE**

Manufactured in our state-of-the-art facilities to the principles of GMP, Klercide Neutral Detergent is safe for use in aseptic core manufacturing areas.



# Klercide™ Neutral Detergent Unit Dose Concentrate



A low foaming, non-ionic detergent solution

## Why use Klercide Neutral Detergent Unit Dose Concentrate?

### PRECISION

Klercide Neutral Detergent is supplied as a liquid concentrate in capped bottles to allow the correct dosing every time. The concentrated liquid has been designed to easily mix, providing a cleaning solution suitable to clean in high-grade cleanrooms.

### SIMPLE

The single unit dose liquid reduces handling and ensures accurate dilution with no filtration required. This results in a cost effective, simple solution to a wide range of cleanroom contamination control challenges.

### DESIGNED FOR CLEANROOMS

Klercide Neutral Detergent is manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These conditions ensure Klercide Neutral Detergent is manufactured to the highest quality and is compliant for use in critical areas.

Klercide Neutral Detergent Unit Dose Concentrate is:

- ▲ Packed and double bagged in ISO class 7 cleanroom
- ▲ Supported with Certificates of Analysis and Sterilization with every batch
- ▲ Supported with full documentation



Klercide Neutral Detergent Unit Dose Concentrate  
76 ml capped

**Precise.  
Simple.  
Efficient.**

### APPLICATION AREA

For use on floors, walls and hard surfaces.

### INSTRUCTIONS FOR USE

Empty all contents into 1 gallon of water of suitable quality for the environment in which it is to be used. Mix gently and when entirely dissolved use immediately.

Use standard procedure to wet surfaces entirely.

### INGREDIENTS

Fatty alcohol ethoxylate water for injection.

### PACKAGING INFORMATION

Klercide Neutral Detergent Unit Dose Concentrate

UNIT SIZE	UNIT OF SALE
76 ml capped	50 x 76 ml

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

### SHELF LIFE

Unopened, this product has a three year shelf life.

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

### WORLDWIDE HEADQUARTERS

1 Ecolab Place  
St. Paul, MN 55102  
USA

[www.ecolab.com/lifesciences](http://www.ecolab.com/lifesciences)  
[www.ecolab.com/pharmaceutical](http://www.ecolab.com/pharmaceutical)