



## IMPREGNATED WIPES FOR EFFECTIVE REMOVAL OF SURFACE CONTAMINATION



### Klerwipe™ Neutral Detergent Pouch Wipes

Non-woven, low particulate wipes impregnated with a neutral, non-ionic detergent in a convenient resealable pouch

**Klerwipe Neutral Detergent pouch wipes are convenient, non-ionic detergent wipes, ideal for routine use of hard surface cleaning and to remove contamination in critical areas.**

Effective on oily or sticky spillages or following general maintenance, and suitable for surface cleaning in pharmaceutical manufacturing areas.

These highly absorbent wipes allow rapid cleaning of surfaces, equipment and production components, leaving surfaces free of soiling to enable effective subsequent sanitization.

Klerwipe Neutral Detergent pouch wipes are:

- ▲ Hydroentangled polyester/cellulose
- ▲ Impregnated with Klercide Neutral Detergent
- ▲ Non-woven and binder free
- ▲ Double bagged in a resealable pouch
- ▲ Basis weight: 58-78 gsm

#### PRODUCT QUALITY ASSURANCE

Manufactured in our state-of-the-art facilities to the principles of cGMP, Klerwipe Neutral Detergent pouch wipes are safe for use in critical aseptic core manufacturing areas.



# Klerwipe™ Neutral Detergent Pouch Wipes



Non-woven, low particulate wipes impregnated with a neutral, non-ionic detergent in a convenient resealable pouch

## Why use Klerwipe Neutral Detergent Pouch Wipes?

### PATIENT SAFETY

Ecolab Klerwipes are absorbent and remain wet throughout use, ensuring a consistent application of Neutral Detergent to the entire surface area. The use of wiping can increase cleaning consistency and the quality of the final product, thereby protecting patient safety.

### EASE OF USE

Ready-to-use, quick and easy to apply, Klerwipe Neutral Detergent wipes are impregnated with Klercide Neutral Detergent, a low-foaming, non-ionic solution with neutral pH. These non-corrosive wipes are powerful and effective in a wide range of cleanroom applications.

Klerwipe Neutral Detergent wipes come in a resealable pouch and are double bagged. The bagging utilizes Klear-tear technology to ensure the bagging can easily be removed during material transfer when wearing gloves. This design also ensures that opening and resealing the pouch does not generate particles, which helps reduce risk in critical areas. In addition, each resealable pouch contains 15 wipes, which helps reduce waste.

### PRODUCT QUALITY

Klerwipes are manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klerwipe Neutral Detergent wipes are manufactured to the highest quality and are compliant for use in your critical areas.

The Klerwipes are:

- ▲ Impregnated with Klercide Neutral Detergent
- ▲ Made with Water for Injection (WFI)
- ▲ Packed and double bagged in an ISO Class 7 cleanroom
- ▲ Supported with Certificates of Analysis and Sterilization with every batch
- ▲ Supported with full documentation



Klerwipe Neutral Detergent  
200 x 200 mm

### WORLDWIDE HEADQUARTERS

1 Ecolab Place  
St. Paul, MN 55102  
USA

[www.ecolab.com/lifesciences](http://www.ecolab.com/lifesciences)  
[www.ecolab.com/pharmaceutical](http://www.ecolab.com/pharmaceutical)

PHARMACEUTICAL

## Safety. Compliance. Efficiency.

### APPLICATION AREA

For use on hard surfaces in all cleanroom areas.

### INSTRUCTIONS FOR USE

For optimum results, ensure complete wetting of surfaces. Re-seal package when not in use. Dispose of package only when emptied and closed.

### INGREDIENTS

Fatty alcohol ethoxylate, tetrahydroxypropyl ethylenediamine and Water for Injection.

### PACKAGING INFORMATION

Neutral Detergent Pouch Wipes

SUBSTRATE	UNIT SIZE	UNIT OF SALE
45% Polyester/ 55% Cellulose	200 x 200 mm wipe	20 packs x 15 wipes

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

### SHELF LIFE

Unopened, this product has a one year shelf life.

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

