

PHARMACEUTICAL



THE QUICK AND EASY ALTERNATIVE TO MOP AND BUCKET SYSTEMS



Klerwipe™ 70|30 IPA Mop Wipes

Non-woven, low-particulate wipes impregnated with 70% v/v Isopropyl alcohol blended with Water for Injection

Ecolab's Klerwipe 70|30 IPA mop wipes are a quick and convenient method for cleaning large surfaces. They are designed for use with the Klercide autoclavable mop wipe frames for optimum surface coverage and remove the need to use space-consuming buckets or trolley systems.

Klerwipe 70|30 IPA Mop Wipes are:

- ▲ Hydroentangled polyester/cellulose
- ▲ Impregnated with Klercide 70|30 IPA
- ▲ Non-woven and binder free
- ▲ Basis weight: 58-78gsm

PRODUCT QUALITY ASSURANCE

Manufactured in our state-of-the-art facilities to the principles of cGMP, Klerwipe 70|30 IPA mop wipes are safe for use even in critical aseptic core manufacturing areas.



ECOLAB®

Everywhere It Matters.™

Klerwipe™ 70|30 IPA Mop Wipes



Non-woven, low-particulate wipes impregnated with 70% v/v Isopropyl alcohol

Why use Klerwipe 70|30 IPA Mop Wipes?

EASE OF USE

Ready-to-use, quick and easy to apply, Klerwipe 70|30 IPA Mop Wipes are impregnated with Klercide 70|30 IPA. These wipes are powerful and effective in a wide range of cleanroom applications.

Klerwipe 70|30 IPA Mop Wipes are designed to fit easily onto the Klercide mop wipe frames. Compared to bucket systems, they are:

- ▲ Faster to use
- ▲ Easier to manage
- ▲ Able to eliminate the potential risk of a dilution error when mixing a concentrate

ENVIRONMENTALLY FRIENDLY

The use of Klerwipe 70|30 IPA Mop Wipes, rather than bucket and reservoir systems helps to:

- ▲ Reduce wasted product
- ▲ Eliminate unnecessary mixing stages
- ▲ Reduce disposal risk of chemistry entering the drainage system

PRODUCT QUALITY

Klerwipes are manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klerwipe 70|30 IPA wipes are manufactured to the highest quality and are compliant for use in your critical areas.

The Klerwipes are:

- ▲ Impregnated with Klercide 70|30 IPA
- ▲ Made with Water for Injection (WFI)
- ▲ Packed and bagged in an ISO Class 7 cleanroom
- ▲ Supported with Certificates of Analysis and Certificates of Irradiation with every batch
- ▲ Supported with full documentation



Klerwipe 70|30 IPA
420 x 250 mm

PHARMACEUTICAL

Safety. Compliance. Efficiency.

APPLICATION AREA

For use on floors, walls and other large, hard surfaces.

INSTRUCTIONS FOR USE

For use with the Ecolab Klercide mop wipe frames. Pull back the upper layer of the peel-pouch pack and place wipe frame head on the mop wipes, wrap the wipe around the head and secure either by pushing firmly into plugs, or with stainless steel clips.

INGREDIENTS

70% v/v Isopropyl alcohol blended with Water for Injection.

PACKAGING INFORMATION

70|30 IPA Mop Wipes

SUBSTRATE	UNIT SIZE	UNIT OF SALE
45% Polyester/ 55% Cellulose	420 x 250 mm Wipe	10 packs x 20 wipes

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

SHELF LIFE

Unopened, this product has a two year shelf life.

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences
www.ecolab.com/pharmaceutical