



Protocol System

Engineered for Cleanroom Environments

Superior Product for Unique Needs

Protocol takes Ecolab's superior cleaning technology and refines it for a cleanroom environment. Adhering to strict cleanroom process control, Protocol delivers the lowest practical particle count while providing exceptional, cost-saving results through consistent finished fabric quality.

Protocol 1000

High-Active Builder

A highly active, ultra pure alkaline concentrate developed to provide extra soil removal and suspension on heavy soil cleanroom textiles. Protocol 1000 is used simultaneously with Protocol 3000 in the break operation of the heavy soil wash formula.

Protocol 3000

Concentrated Neutral Detergent

A super concentrated, ultra pure detergent designed to penetrate, emulsify and lift off all types of cleanroom soils. It disperses rapidly, instantly wetting and penetrating the fabric for quick, deep soil removal.

Protocol 4000

Ultraconcentrated Sour

Safe on all cleanroom fabrics, this concentrated neutralizing agent is formulated with a blend of mild acids for maximum neutralization capacity.

As the global laundry care expert, only Ecolab offers you comprehensive solutions that help optimize your entire operation. Partner with us to get the proactive service you need to gain the competitive advantage you want. Our outstanding territory managers and service specialists offer individual facility consultations. We assess your unique needs and customize wash formulas to deliver exceptional results.



Talk to your Ecolab
Territory Manager or contact us
at 1.800.553.8683

Meeting Needs, Big and Small

Ecolab's Textile Care Division provides products, systems and services that meet the complex needs of large commercial laundry operations. Ecolab recognizes that each customer differs in size, technology and service requirements and delivers a total approach to laundry problems with a wide and varied offering that is unmatched in the industry.

Value-Added Services

Ecolab focuses sales, research and development resources in ways that are exclusively targeted to give you a competitive edge. Ecolab associates understand the challenges commercial laundriers face. Not only do our products and technology help produce greater results, but our industry expertise can also add great value to your organization.

Experienced Laundry Specialists

Ecolab's Textile Care territory managers are veteran problem-solvers, trained specifically for the complexities of commercial laundries. They know tough commercial soils, heavy-duty equipment and the formulas needed for soil removal, energy savings and washaisle efficiency. Most importantly, they understand that every laundry is different and are trained to custom-design wash programs accordingly.

Superior Solutions, Year After Year

Results are everything when it comes to laundry operations. Labor, equipment, chemicals and utilities all impact the bottom line. Partner with Ecolab to get the results you need, so you can focus on growing your business.

